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
Achieves simultaneous multipoint photoactivation

The new DMD module enables photoactivation and photoconversion of a user-specified pattern and position(s), whereas the conventional FRAP unit only enables photoactivation of a single, manually-positioned spot. The DMD illumination shape, size, position and number can be freely customized using the NIS-Elements software. This capability allows researchers to optically mark a subset of cells or protein populations within a single cell or multiple cells to track their behavior. The DMD module is also optimally suited for optogenetics experiments in which highly customized ROIs can be

used to optically induce functional changes in subsets of cells or protein populations. The DMD module can be used with either laser illumination or less phototoxic LED illumination.

A mouse embryonic fibroblast co-expressing mCherry-tagged lamin A (red) and photo-activatable GFP-tagged lamin A was photo-converted (green) in the lower right region using the DMD module and 405 nm LED light. Time-lapse images were captured using the epi-fluorescence illuminator. By photoactivating a sub-population of the lamin proteins, one can observe their dynamics and subunit-exchange behavior.

Image courtesy of Drs. Takeshi Shimi and Bob Goldman, Northwestern University Medical School

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
Fully automated TIRF adjustment and observation is now possible

The ideal incident angle and focus of the laser for TIRF observation vary depending on specimen and observation conditions. Adjusting the incident angle and focus for achieving TIRF requires skill and experience. The new H-TIRF module automatically adjusts the focus and incident angle of the laser for TIRF observation by monitoring the reflection beam. This automatic laser focus adjustment and incident angle adjustment is carried out by the auto-alignment function in NIS-Elements software. Incident angles and penetration depths of the evanescent fields can be saved and reproduced for subsequent experiments to ensure consistent imaging results.

The H-TIRF module is configured with a gradation neutral density (ND) filter that can be moved into the light path to achieve an even field of TIRF illumination.

An in vitro preparation of fluorescently-labeled microtubules (tetramethylrhodamine and Alexa 647) and tubulin binding proteins (Alexa 488) was imaged in three different wavelengths using the H-TIRF illuminator and the gradation ND filter. Incident angles can be automatically adjusted for multiple wavelengths.

Image courtesy of Melissa Hendershott and Dr. Ron Vale, University of California, San Francisco

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Camera Flash 4 (Hamamatsu)

Quantum efficiency	Over 70 % at 600 nm and 50 % at 750 nm
Imaging device	Scientific CMOS Sensor FL-400
Effective no. of pixels	2048 (H)×2048 (V)
Cell size	6.5 μ m×6.5 μ m
Effective area	13.312 mm×13.312 mm
Full well capacity	30 000 electrons (typ.)
Readout speed	100 frames/s (Full resolution, Camera Link)

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